

# Johnson & Johnson Vitros DT60 II Chemistry System

Manufacturer: KNAUER

SKU#: 8195-10-0020

- Extensive Test Menu
- Quick Results
- Dry-Slide Convenience
- Small Sample Size
- Efficient Operation



### **Product Description**

Johnson & Johnson's Vitros DT60 II chemistry system efficiently offers accurate diagnostic test results. This cost-effective model is an easy-to-use bench top model for strict point-of-care and back-up testing requirements. Laboratory needs including credible test results, reduced overall costs and increased patient satisfaction are conveniently met by Vitros DT60 units. Users can efficiently operate DT60 II systems with the once-a-day control adjustment and twice-a-year calibration feature to save valuable time. Operators can ensure the reliability and productivity of laboratory procedures with this chemistry system. Extensive Test MenuThe Vitros DT60 chemistry system is capable of performing a broad range of chemical applications including basic chemistries, electrolytes, enzymes, full lipid profile and derived tests. The system's random access testing capability allows users to run chemistries of their choice in any order. Quick Results Users can take advantage of the quick sample processing of Vitros DT60 II systems, obtaining results and hard copy printouts with patient ID in as fast as two minutes. High throughput rates also allow the delivery of up to 10 tests in five minutes. Approximately 100 tests can be processed within an hour. For best results, the use of small sample sizes per test is encouraged and facilitates testing hard-todraw geriatric and pediatric patients. Fingerstick or venipuncture is required in obtaining 10 microliters in every cycle.Dry-Slide Convenience VITROS dry-slide technology minimizes the effects of common interferences. This technology makes sure that results are precise while increasing the system's cost and convenience benefits. Users can avoid mixing wet reagents, wasting samples as well as using special plumbing with this model's integrated dry-slide technology. Call us today to learn more about the Johnson & Johnson Vitros DT60 II Chemistry System!

### Johnson & Johnson Vitros DT60 II Chemistry System Specification

## Johnson & Johnson Vitros DT60 II Chemistry System Specifications

DT60 II Analyzer

 Width
 47.6 cm (18.75 in)

 Depth
 34.9 cm (13.75 in)

 Height
 17.1 cm (6.75 in)

Weight Approximately 8.6 kg (19 lbs.)

DTSC II Module

 Width
 34.3 cm (13.5 in)

 Depth
 35.4 cm (13.9 in)

 Height
 16.5 cm (6.5 in)

 Weight
 7.7 kg (17 lbs)

DTE II Module

1

Width Depth Height Weight Voltage for all components 14.6 cm (5.75 in) 35.4 cm (13.9 in) 16.5 cm (6.5 in) Approximately 2.7 kg (6 lbs.)

110 - 240 VAC (50-60 Hz) RS-232 port for computer interface

#### WHY BUY FROM GMI?

Your trusted source for Used & New lab equipments, Parts, Factory Clearance instrument service since 1985.







Shop a comprehensive selection of brand new lab instruments from trusted brand purchases with our Trade In & Trade Up offer training and installation options as names.

Enjoy big savings on New Instrument Program.

Experience customer-driven services. We well as multi-year service agreements.

GMI is a ISO 9001:2008 certified company - your proof that your instrument purchase from GMI passes through rigorous refurbishing, recalibration, and recertification process, arriving at your laboratory in fully working order and defect free. We are not an auction company, nor are we an instrument broker... getting in the middle of 'deals'. All items come complete with the final QC date, engineer's ID# and a copy of test results / QC data, if applicable. Additionally, all items are thoroughly cleaned an professionally packaged / crated with care to ensure they arrive safely in your lab ready for professional installation.

